CEMENTO Y C~Plug®

Cemento óseo y restrictor reabsorbible
Bone cement and resorbable restrictor

Cemento óseo Bone cement

Libertad para elegir

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Freedom to choose

In 1958 the first implantation of a total hip endoprosthesis using a PMMA-cement took place. Today, nearly 50 years later, the use of bone cement is a common procedure in total hip and knee arthroplasty. It is without question, that the cements of “the first hour” have a backup of performance, excellence and proven clinical success.

These days, more than 40 different bone cements are commercially available. The surgeon can choose between low, mid and high viscosity cements. The products are available with and without antibiotics. All the cements are, of course, within the specification according to ISO 5833 criteria and all are in clinical use.

There must be more to find the right cement

Beside all the technical and chemical properties of different bone cements, the working behaviour of a bone cement is an important issue to distinguish a good cement.

The current composition of the different versions of the line has taken into account the long acquired experience of major cements available. Offers a high quality without losing the well known properties of proven components.

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The composition in particular

Buthylmethacrylate (BMA)

Buthylmethacrylate (BMA) is used as a copolymer. It is well known, that BMA causes a reduction of bone cement setting temperature by lowering the exothermic polymerization reaction, while generating fewer vapours. It is also founded that cements with BMA will transfer the stress quickly and distribute it over a larger area of endosteal bone surface. The so improved elasticity enhances long term results.

Barium sulphate

Radiopacity is needed to visualize bone cement on X-ray examination. Commercially two radiopaque elements are in use: Bariumsulfate and Zirconiumdioxide.

The choice of barium sulphate as a radiopacifier is also a rational choice. Literature and studies have shown, that barium sulphate causes less damage to bearing surfaces and lower wear rates while reducing particulate debris.

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
Low setting temperature

Polymerization is an exothermic chemical reaction. During its hardening phase the temperature of the cement increases.

The thermal effect of the polymerization process enhances necrotic effects. For that reason the setting temperature of bone cements should be within the ISO 5833 (2002) limit of 90° C degrees.

The peak exothermic temperature for the Cemento line lies between 67° C (Genta Cemento 1) and 79° C (Genta Cemento 3).

Both values are clearly lower than the ISO limit and also lie under the setting temperatures of competitive products.

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
Deep intrusion for a long lasting result

Polymerization is an exothermic chemical reaction. During its hardening phase the temperature of the cement increases.

The thermal effect of the polymerization process enhances necrotic effects. For that reason the setting temperature of bone cements should be within the ISO 5833 (2002) limit of 90° C degrees.

The peak exothermic temperature for the Cemento line lies between 67° C (Genta Cemento 1) and 79°C (Genta Cemento 3).

Both values are clearly lower than the ISO limit and also lie under the setting temperatures of competitive products.

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
Two Viscosities for your preference


After mixing the powder and the liquid components, a doughy mass is formed by swelling and dissolution. Because of the initiated radical polymerization the viscosity increases continuously until the complete curing of the cement. Viscosity is an important handling property for the surgeon.

Low, medium and high viscosity cements can be distinguished.

Cemento 3 is a low viscosity bone cement which is initially more liquid and therefore very well suited for vacuum mixing without the need of pre-cooling.

The high viscosity cements Cemento 1 are initially ductile. Handmixing and knotting are recommended.

With Cemento line a variety of cements is available which meets the requirements of various surgical procedures and various surgeon's preferences and techniques.

Cemento 1 is recommended for: Total Knee, Hybrid Knee, Uni Knee, Shoulder Stem, and Glenoid.

Cemento 3 is recommended for: Hip, Hip-Resurfacing, Total Knee, Hybrid Knee, Uni Knee, Elbow, Shoulder Stem, and Glenoid.

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
Similar properties within the whole product range


As mentioned earlier, it is not surprising that all 4 cements comply with the ISO 5833 (2002). One goal of the development of Cemento line was to reach a high, persistent quality with comparable mechanical properties of all products.

It is remarkable that even with different handling behavior all values are within one range, without significant differences.

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Antibiotics in bone cements


For powdered antibiotics the critical value of admixture is about 3-4 %. If more antibiotic powder is added, the density of the bone cement matrix and the compressive and bending strength could be influenced. The amount of an antibiotic is a compromise between the antibiotic concentration at the target site and the influence on the mechanical properties of the bone cement. Therefore the amount of antibiotic added to Genta Cement 1 and Genta Cement 3 is 1,2%. Consequently, the effect on the mechanical strength is neglectable.

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Release of gentamicin

In Europe, the most frequently used bone-cements are loaded with Gentamicin. Gentamicin covers 75% of the current bacteria.

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
The release of Gentamicin from bone cement is related to the cement on the one hand, and to the type and dosage of the antibiotic on the other. The admixture of antibiotics to the cement matrix leads to local release of the antibiotic at the site of implantation. The initially high antibiotic concentration prevents the adherence of bacteria to the cement surface leading to colonization and infections.

The release rates of Genta Cement correspond to those obtained from other Gentamicin loaded bone cements with high initial and then steadily decreasing concentrations.

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BONE CEMENT

GIVES THE SURGEON THE FREEDOM TO CHOOSE

- There is more than one cementing technique

PROVEN AND WELL ESTABLISHED PRODUCT

- 10 years of clinical experience
- More than 200.000 cases without complaints

EXCELLENT HANDLING CHARACTERISTICS

- The right products for hand - or vacuum mixing, without pre-cooling
- Efficient working time and short waiting time
- Powder in a bottle

THE RIGHT CEMENT FOR THE CHOSEN TECHNIQUE

- PMMA Bone Cement with two different viscosities meets the requirements

CONSISTENT QUALITY

- High degree of reproducibility

Cemento óseo
Bone cement

Ref. F9000001	Cemento 1
Ref. F9000003	Cemento 3
Ref. F9000011	Cemento Genta 1
Ref. F9000013	Cemento Genta 3

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A modern and safe plug should fulfil some requirements

- Biodegradability
- Biocompatibility
- Easy Instrumentation
- High Pressure Resistance

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But C~Plug® is showing more:

C~Plug® is a biodegradable and biocompatible bone cement restrictor out of 50% gelatine, 30% glycerol and 20% water that is available in 6 sizes (from 8mm to 18mm).

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
1. Accurate positioning

The insertion hole of C~Plug® is slightly conical, the insertion instrument straight. The result is less suction between the instrument and the plug when the instrument is withdrawn – the plug is not pulled backwards.

The extra depth of the insertion hole ensures a stable connection with the applicator during the setting, capsizing of the restrictor is avoided.

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
2. Increasing the contact surface

The basic concept of C~Plug® is a cylindrical plug with three rings seated on the base. The specially designed rings lead to an increasing of the contact surface during insertion. This type of curved rings has a high degree of resistance to displacement, a great stability and supports the bone cement penetration perfectly.

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
3. Reducing the distal pressure

With some current resorbable restrictors, a very high distal pressure is generated.

Due to the outstanding PPRS® (Peak Pressure Release System), the distal pressure with C~Plug® is reduced to 300 millibars. A special shaped air vent allows the air to pass. It also helps to balance the pressure on both sides of the device: the C~Plug® stays exactly where initially placed.

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4. All-in-one-application System

The innovative design of the instrument set allows to measure the femoral canal and to insert the C~Plug® with same device.

The insertion rod is scaled, allowing for pre-insertion depth calculation and the plug carriers provide exact diameter identification.

A longer screw thread and a solid inserter make the system a stable and long term usable device.

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The specific concept C-Plug[®] with its elastic and mechanical properties results in:

- Homogenous cement mantle
- Deep penetration of the cement
- Efficient protection of the medullary canal
- Reduced distal pressure
- Total resorption within some days
- No plug removal for revision surgery


Restrictor de cemento óseo reabsorbible C-Plug [®] C-Plug [®] resorbable bone cement restrictor			
Ref. F9000108	8 mm	1 pieza	1 pcs
Ref. F9000110	10 mm	1 pieza	1 pcs
Ref. F9000112	12 mm	1 pieza	1 pcs
Ref. F9000114	14 mm	1 pieza	1 pcs
Ref. F9000116	16 mm	1 pieza	1 pcs
Ref. F9000118	18 mm	1 pieza	1 pcs
Ref. F9000120	Set instrumental	1 varilla de inserción y 6 calibradores	1 inserter & 6 sizes in steri-box

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IDICENCA01 / 2011 ©