

TiN - Titanium Nitride



*Gold Standard
the strong coating...*

implantcast 



e

ic coating that has many advantages

increased wettability

extreme hardness

cemented and cementless components

high wear-resistancy

reduced ion release

excellent adhesive strength

minimized polyethylene wear

reliable coating

more than 20 years of experience

more than 300.000 implants without flaking



Reliable coating


For more than 20 years titanium nitride (TiN)-coatings have been applied to prosthesis to prevent implant allergies and to reduce wear in total knee arthroplasty. TiN-coated implants have been established successfully in the market.^{1, 2, 3, 20} Ideal hygienic conditions in the manufacturing environment as well as the correct coating type and process management ensure constant product safety. Implantcast GmbH cooperates with DOT GmbH, Rostock (Germany) a company specialized on coatings with more than 20 years of experience in coating medical implants.⁴ Every implant coating is 100% inspected ensuring the highest possible quality standards. More than 300.000 TiN-coated cemented and cementless implants have been clinically successfully applied to date.

Base material

Base for a coating with titanium nitride is a finished implant with a highly polished articulating surface. The implant components of the ACS®- system are manufactured from implavit®, a cobalt-chromium-molybdenum (CoCrMo)-alloy according to ISO 5832-4. In the process of manufacturing carbides (bondings of metal and carbon) develop, which are characterized by extreme hardness and these carbides partly protrude from the surface of the finished implant at a microscopic level. These protruding carbides cause wear of the softer articulating bearing surface - in this case polyethylene.

element	cut-off grade % (weight per- cent)
Chromium	26,5-30
Molybdenum	4,5-7
Nickel	max. 1,0
Iron	max. 1,0
Carbon	max. 0,35
Manganese	max. 1,0
Silicium	max. 1,0
Cobalt	rest

CoCrMo-alloy
chemical composition
acc. to ISO 5832-4



Structure of a highly polished non-coated component


cemented CoCrMo-component coated with TiN

TiN covers carbides

A limiting factor for the life time of a knee replacement is its cyclic and load dependent wear of the softer polyethylene after implantation which produces a significant release of wear debris. The carbides that protrude at the implant surface, which are mainly responsible for wear of the polyethylene, are covered by the much harder TiN-coating. This results in severely reduced wear of the softer articulating partner.⁵

Technology

For manufacturing of the ceramic TiN-coating a specific arc vaporization technique (PVD-coating, physical vapour deposition) is applied. Nitrogen supply is added to the implants as they are coated in an evacuated vacuum chamber in the vapor phase. The computer-controlled process leads to a high reproducibility and coating uniformity. The process anchors the coating safe in several layers of atoms to the implant surface. Therefore only the implant surface is modified. The material properties of the base material (CoCrMo-alloy) as well as its biomechanical functionality remain unaltered.



highly polished non-coated component versus coated component

cementless CoCrMo-component coated with TiN


TiN - key advantages

Reduced wear

Wear of the polyethylene is in the longterm one of the main causes for failure of knee endoprosthesis. Aseptic implant loosening can result due to the process of abrasion and wear and consequently revision surgery becomes inevitable.

The TiN-coating reduces wear to a minimum. In in-vitro wear tests on a knee simulator according to ISO 14243 standard the high wear-resistancy of the hard TiN-coating was proven over CoCrMo.¹ The wear rate with TiN-coated components was 38% of the wear rate of the non-coated CoCrMo-components.⁶

This demonstrates that the TiN coating has superior bearing qualities to CoCrMo.



Wear rates of the ACS® knee system after 5 Mio cycles⁶

In wear tests the ceramic TiN-coating demonstrates a higher resistancy against scratching. Particles of bone cement are temporarily tolerated in the tribological pairing because extremely hard foreign particles merely generate striae. Thus the potential risk of third body wear is minimized. On the contrary at a non-coated implant third bodies can induce small scratches in the surface, which lead to increased wear of the polyethylene.



Surfaces after Pin-on-Disc-Test with bone cement, TiN (left), non-coated titanium (right)



Minimized ion release

All metallic implant components release ions to their environment over time. In some patients such ions can elicit allergic reactions. Nickel, cobalt and chromium, which are elements of the base material of the articulating implant components, are considered the most frequently allergy eliciting metals.¹⁹


The TiN-coating is biocompatible and acts like a barrier; the potential release of allergy eliciting ions of the base material is reduced to a minimum.⁷ The ceramic coating itself is inert in the human body.

In in-vitro tests TiN-coated CoCrMo was suspended and after a defined time of exposure the ion concentration of the solution was analyzed with regard to Co-, Cr- and Ni-ions. This was compared to the data from non-coated CoCrMo. There was a significant difference in metal ion release in favour of TiN coating.

The ion release is so minimal in the TiN coated components that it is below the limit of detectability.⁸

Also in clinical practice there have never been any evidence of allergic reactions with implants that have been TiN-coated showing an intact surface.³

Therefore the TiN-coating on implant components is especially suitable for patients with sensitivity to nickel, chromium or cobalt.^{3, 8}



Metal ion release for TiN-coated and non-coated CoCrMo⁸

TiN and its properties

Adhesive strength

The adhesive strength is one of the most important parameters for the quality of a coating coupled with how well it adheres to the base material. The coating is classified into adhesive strength categories (HF), in which HF1 means an excellent adhesive strength and HF6 an insufficient adhesive strength that produce flaking of the coating from the base material.

TiN-coating onto a CoCrMo base material has an adhesive strength of HF1.⁹


In in-vitro tests there is no evidence of the TiN-coating flaking (delaminating) from the base material.⁹ The adhesive strength of the TiN-coating is technologically ruled.^{1, 2} Optimal hygienic conditions in the manufacturing environment as well as the correct type of coating and process management ensure a constant product safety.

Previous studies have proven that the adhesive strength of TiN-coatings on CoCrMo meets the top requirements and it is rated higher than the adhesive strength of a zirconium oxide diffusion layer (Oxinium) on a ZrNb-alloy (HF1 vs. HF3).¹⁰

In the past it has been reported that some clinical partial failure of the coating occurred on femoral heads.^{11,12} The so-called „eggshell effect“ was held responsible for that. This means that the coating should flake off the soft base material when interacting with third bodies (particles) in the presence of high punctual loads (bearing pressure).

However on closer examination of these failed cases in specific regard to coating adhesion strength, it was found that they were manufactured by an inferior manufacturing process that was widely adopted at that time. As a consequence they exhibited an insufficient adhesive strength and partial failure occurred.^{1, 8} The manufacturing process pioneered by DOT GmbH has never had a reported case of delamination.

DOT GmbH have more than 15 years without any report on adhesion failures. The adhesive strength of the coating is so strong that even particles of bone cement are tolerated without problems in the tribological pairing. Extreme hard foreign particles merely generate striae on the surface, which however do not result in delamination of the coating.⁴




Determination of the adhesive strength: HF1 (left) and “eggshell effect” HF6 (right)

Rockwell-Test according to DIN 50103:

A diamond pin is pushed in the base material with a defined load. Afterwards the coating of the boundary is optically classified in the adhesive strength categories HF1-HF6 according to scratches and flaking of the coating.

Surface hardness

The hardness of the surface plays an important factor in the resistance against wear. The surface hardness of the TiN-coating is 23,5 GPa. This is 5 times higher than that of a non-coated CoCrMo-component.





Hardness of different surfaces^{1,13}

Wettability

When a water droplet is placed on given material the contact angle shows the wettability of a surface which directly correlates to the gliding properties and surface lubrication of a material. The smaller the contact angle between surface and fluid, the better is the wettability and the hydrophilic property of the surface.



A drop of fluid on a TiN-coated implant covers a larger surface area compared to a drop of fluid on a non-coated CoCrMo-component. This increased wettability of a TiN-coated surface when in-vivo surrounded by synovial fluid reduces the friction between femoral component and PE-insert and therefore reduces the wear of the PE-insert.



Wettability of a TiN-coated and non-coated CoCr-component with bovine serum⁹

Adhesive strength of bone cement on TiN

In-vitro tests have proved that TiN coating has a superior adhesion strength of bone cement over CoCrMo. The tests demonstrated a significant increase of one third of the adhesive strength of bone cement on TiN-coated surfaces compared to non-coated surfaces.²¹




Adhesive strength of bone cement on different surfaces²¹

TiN in comparison...

...to other implant surfaces


There are alternative modified surfaces of endoprosthesis in the market, which also have been designed to overcome abrasive wear as well as minimized ion release. These include amongst others multi-layer ceramic coatings (AS-coating, Aesculap) and an oxidized surface of the implant (Oxinium, Smith and Nephew).

TiN




- monolayer coating
- ceramic
- coating thickness app. 5,5µm
- surface hardness 23,5GPa
- coating colour golden

AS¹³



- multilayer coating (7 layers)
- ceramic
- coating thickness total app. 4,5µm
- surface hardness 25GPa
- coating colour golden

Oxinium^{13, 14, 15}



- zirconium-niobium-metal alloy with oxidized surface
- ceramic
- no coating
- thickness oxidized surface app. 5µm
- surface hardness 12GPa
- colour black

The mean coating thickness of the TiN-coating (app. 5,5µm) is greater than the thickness of the multilayer AS-coating.¹³ This increased coating thickness ensures increased safety against wear in case of emerging third body wear by particles of bone cement for instance.

In regard to Oxinium - if the oxidized surface fails by wear (e.g. third body wear by particles of cement), there will be just soft zirconium metal as articulating partner. Zirconium is approximately as soft as pure titanium and thus not suitable as a bearing surface as pure material.¹⁰

The multilayer architecture of the AS-surface requires demanding process management. Thus its manufacturing process is more prone to small technological failures, which eventually could have an influence on the quality of the coating, compared to the established TiN-coating. An increased adhesive strength by the multilayer structure could not be proven in a publication to date.

The surface hardness of the TiN- and AS-coating is comparable, whereas Oxinium exhibits 50% of the surface hardness of TiN and AS.¹³ This makes the Oxinium surface potentially more sensitive to scratch formation and wear.

All 3 surfaces minimize the ion release of the relevant metals Co, Cr, Ni in in-vitro tests to or below the limit of detectability.^{16, 13, 17, 18} Nevertheless it needs to be noted that the AS-coating even contains chromium as chromium nitride, thus a potentially allergenic metal, in its coating architecture.

Furthermore all 3 surfaces reduce PE-wear compared to a non-coated CoCrMo-component.^{6, 13, 14, 15, 16, 17} But a direct comparison of wear rates of different manufacturers needs to be considered critical due to different implant designs, implant materials and testing conditions.

References

- 1 Bader R et al. (2008) Alternative material and solutions in total knee arthroplasty for patients with metal allergy. *Orthopäde* 37: 136-142
- 2 Thomas P et al. (2008) Orthopedic surgical implants and allergies. *Orthopäde* 37: 75-88
- 3 Baumann A. (2001) Keramische Beschichtungen in der KTEP Standardlösung für Allergiker. *JATROS Orthopädie & Rheumatologie* 6: 16-17
- 4 Baumann A. (2002) DOT GmbH Technischer Report. Die Haftfestigkeit von Titannitrid (TiN) - Schichten auf Implantatlegierungen. Data on file
- 5 Jones V.C. et al. New materials for Mobile Bearing Knee Prosthesis - Titanium Nitride Counterface Coatings for Reduction in Polyethylene Wear. In: Hamelynck K.J., Stiehl L.B. *LCS-Mobile-Bearing Knee Arthroplasty 25 Years of Worldwide Experience* ch21. Springer. New York 2000
- 6 Test report A219/05.1 IMA Dresden. Data on file
- 7 Wisbey et al. (1987) Application of PVD TiN coating to Co-Cr-Mo based surgical implants. *Biomaterials*, 11
- 8 Prof. Thomas LMU München Final Report Effect of a TiNbN or TiN surface coating on cobaltchromium-molybdenum and stainless steel test specimens regarding the release of nickel, chromium and cobalt: evaluation via eluate analysis and in-vitro cytokine release from peripheral human blood cells, Data on file
- 9 Dempwolf H (2013) Mechanical and tribological performance of ZrN and ZrN/CrN PVD layers compared to established PVD coatings in endoprosthetics, Bachelor thesis Universität Rostock
- 10 Galetz MC et al. (2010) Abrasion resistance of oxidized zirconium in comparison with CoCrMo and titanium nitride coatings for artificial knee joints. *Journal of Biomedical Materials Research B: Applied Biomaterials* Vol. 93B Issue 1
- 11 Raimondi MT, Pietrabissa R (2000) The in-vivo wear performance of prosthetic femoral heads with titanium nitride coating. *Biomaterials* 21: 907-913
- 12 Harmann MK et al. (1997) Case report: Wear analysis of a retrieved hip implant with titanium nitride coating. *Journal of Arthroplasty* Vol. 12 No. 8
- 13 Aesculap Knie Arthroplastik Oxinium versus AS Kniesysteme Prospekt Nr. 039101 0211/1.5/1
- 14 Internetseite <http://www.smith-nephew.com/professional/products/all-products/oxinium/oxinium-material/hardness/>
- 15 Smith and Nephew (2010) Long Term Wear White Paper. 30 year wear performance of LEGION™ Primary knees with VERILAST™ Technology
- 16 Aesculap AS Advanced Surface Broschüre Prospekt Nr. 043701 0312/1/1
- 17 Reich et al. (2010) Preclinical Evaluation of coated knee implants for allergic patients. *Orthopäde* Vol 39: 495-502
- 18 Lützner et al. (2013) Metal hypersensitivity and metal ion levels in patients with coated or uncoated total knee arthroplasty: a randomised controlled study. *International Orthopaedics* 37: 1925-31
- 19 Eben R et al. (2009) Implant allergy register - a first report. *Orthopäde* 38: 557-562
- 20 Mohammed et al. (2014) Medium term outcome of Titanium Nitride, mobile bearing total knee replacement. *Acta Orthop. Belg.* Vol. 80: 269-275.
- 21 Radant D DOT GmbH (2013) Project repoort 321-13 Vergleich der Haftfestigkeit von Knochenzement auf unbeschichteten und TiN beschichteten CoCrMo Oberflächen (Data on file)



implantcast

implantcast GmbH
Lüneburger Schanze 26
D-21614 Buxtehude
Germany
phone: +49 4161 744-0
fax: +49 4161 744-200
e-mail: info@implantcast.de
internet: www.implantcast.de



Your local distributor:

TINOVERE-211114

