

HAND & FOOT SYSTEM

WISE-LOCK MODULAR **HAND & FOOT** SYSTEM



about us

Auxein Medical is an integrated, research based, orthopaedic Implants & Instruments manufacturing company, producing a wide range of quality, affordable generic implants, trusted by healthcare professionals and patients across geographies. It is the Company's constant endeavor to provide a wide basket of generic and innovator products that exceed the highest expectations of our customers in term of quality and safety. The company has world-class manufacturing unit established in India and serves customers in over 60 countries worldwide.

Our Achievements

US FDA 510K APPROVED •

(Locking Plates and Screws)

EC Certificate - Class IIa & IIb •

ISO 13485:2016 •

ISO 9001:2015 •

India FDA Approved •

GMP Certificate •

Orthopaedic Implants

Contents

1.5mm Wise-Lock Modular Hand System	3-6
2.0mm Wise-Lock Modular Hand System	7-11
2.4mm Wise-Lock Modular Hand System	12-16
2.7mm Wise-Lock Modular Foot System	17-21
Instrument & Implant Set	22-30
Universal Instruments	31-35

Adaptable Locking Plates for Various Indications

The modular Wise-Lock Hand & Foot System offers a wide range of anatomical and adaptable plates. All Implants are available in Titanium.

Wise-Lock Hand & Foot Implants are intended for the following indications:

1.5mm Wise-Lock Hand System

- Fractures, osteotomies and arthrodeses of the phalanges and metacarpals
- Replantations and reconstructions of the phalanges and metacarpals

2.0mm Wise-Lock Hand System

- Fractures of the phalanges, metacarpals, wrist bones and distal radius (double-plate technique)
- Osteotomies and arthrodeses of the interphalangeal joints


2.4mm Wise-Lock Hand System

- Fractures of the metacarpals, wrist bones and distal radius (double-plate technique)
- Osteotomies and arthrodeses on the hand

2.7mm Wise-Lock Foot System

- Osteotomies and arthrodeses on the foot






4 Modules, 4 Dimensions:

- Adaptable plates for different indications
- Complete range of implants

1.5mm Wise-Lock
2.0mm Wise-Lock
2.4mm Wise-Lock
2.7mm Wise-Lock

- Broad range of locking and Wise Lock plates to address various fracture patterns
- The use of locking screws allows creating a fixed-angle construct, which is particularly advantageous in osteopenic bone or in fractures near the joint



1.5mm Wise-Lock Plate, Straight

Code	Holes	Length
TI-554.04	4	23mm
TI-554.06	6	36mm



1.5mm Wise-Lock Adaption Plate, Straight

Code	Holes	Length
TI-555.06	6	29mm
TI-555.12	12	59mm



1.5mm Wise-Lock Adaption T-Plate

Code	Holes	Length
TI-556.08.03	Shaft 8, Head 3	44.5mm
TI-556.08.04	Shaft 8, Head 4	44.5mm



1.5mm Wise-Lock Strut Plate

Code	Holes	Length	Type
TI-823.00R	8	23mm	Right
TI-823.00L	8	23mm	Left



1.5mm Wise-Lock Adaption Y-Plate

Code	Holes	Length
TI-558.08	Shaft 8, Head 3	47mm



1.5mm Wise-Lock Condylar Plate

Code	Holes	Length
TI-559.06	Shaft 6, Head 2	36mm
TI-559.08	Shaft 8, Head 2	46mm



1.5mm Wise-Lock Screws, Self Tapping, (Star Head)

Code	Length
TI-1130.06	6mm
TI-1130.07	7mm
TI-1130.08	8mm
TI-1130.09	9mm
TI-1130.10	10mm
TI-1130.11	11mm
TI-1130.12	12mm
TI-1130.13	13mm
TI-1130.14	14mm
TI-1130.15	15mm
TI-1130.16	16mm
TI-1130.18	18mm
TI-1130.20	20mm
TI-1130.22	22mm
TI-1130.24	24mm



1.5mm Cortical Screws, Self Tapping, (Star Head)

Code	Length
TI-1281.06	6mm
TI-1281.07	7mm
TI-1281.08	8mm
TI-1281.09	9mm
TI-1281.10	10mm
TI-1281.11	11mm
TI-1281.12	12mm
TI-1281.13	13mm
TI-1281.14	14mm
TI-1281.15	15mm
TI-1281.16	16mm
TI-1281.18	18mm
TI-1281.20	20mm
TI-1281.22	22mm
TI-1281.24	24mm



1.8mm Buttress Screws, (Star Head)

Code	Length
TI-1103.112	12mm
TI-1103.114	14mm
TI-1103.116	16mm
TI-1103.118	18mm



1135-01

Double Drill Guide 1.5/1.1mm


1135-02

Wise-Lock Drill Sleeve 1.5mm For Drill Bit $\phi 1.1$ mm

1135-03

Bending Pin With Thread For 1.5mm Wise-Lock Plates


1135-04

Drill Bit Quick Coupling $\phi 1.1$ mm x Length 75/61mm, 2-Flute

1135-05

Drill Bit Quick Coupling $\phi 1.5$ mm x length 65/50mm, 2-Flute

1135-06

Countersink For 1.5mm Screws



1135-07

Screwdriver Shaft, T4 for 1.5mm Screws


1135-08

Screw Holding Sleeve for 1.5mm Cortical Screws


1135-09

Screw Holding Sleeve for 1.5mm Wise-Lock Screws


1138-08

Container with Tray and Lid For 1.5mm Implant & Instrument Set



2.0mm Wise-Lock Plate, Straight

Code	Holes	Length
TI-505.204	4	31mm
TI-505.205	5	38mm
TI-505.206	6	45mm
TI-505.207	7	52mm
TI-505.208	8	59mm
TI-505.210	10	73mm



2.0mm Adaption Plate, Straight

Code	Holes	Length
TI-821.212	12	71mm



2.0mm Wise-Lock Adaption Plate, Straight

Code	Holes	Length
TI-555.212	12	81.26mm



2.0mm Wise-Lock Condylar Plate

Code	Holes	Length
TI-559.204	Shaft 4, Head 2	33.6mm
TI-559.205	Shaft 5, Head 2	40.4mm
TI-559.206	Shaft 6, Head 2	47.2mm
TI-559.207	Shaft 7, Head 2	54.0mm



2.0mm Wise-Lock Rotation Correction Plate

Code	Holes	Length
TI-989C.204	Shaft 4, Head 2	34mm
TI-989C.205	Shaft 5, Head 2	40mm



2.0mm Wise-Lock T-Plate

Code	Holes	Length
TI-556.204	Shaft 4, Head 3	32.1mm
TI-556.205	Shaft 5, Head 3	38.9mm
TI-556.206	Shaft 6, Head 3	45.7mm
TI-556.207	Shaft 7, Head 3	52.5mm



2.0mm Wise-Lock T-Adaption Plate

Code	Holes
TI-555.207T	Shaft 7, Head 2



2.0mm Wise-Lock Y-Adaption Plate

Code	Holes
TI-555.207Y	Shaft 7, Head 3



2.0mm Mini Condylar Plate

Code	Holes	Overall Height	Type
TI-515.07L	7	29mm	Left
TI-515.07R	7	29mm	Right



2.0mm Mini H-Plate

Code	Holes	External Dimensions
TI-814.04	4	11x12mm



2.0mm Cortical Screws, Self Tapping, (Star Head)

Code	Length
TI-102.206	6mm
TI-102.207	7mm
TI-102.208	8mm
TI-102.209	9mm
TI-102.210	10mm
TI-102.211	11mm
TI-102.212	12mm
TI-102.213	13mm
TI-102.214	14mm
TI-102.215	15mm
TI-102.216	16mm
TI-102.218	18mm
TI-102.220	20mm
TI-102.222	22mm
TI-102.224	24mm
TI-102.226	26mm
TI-102.228	28mm
TI-102.230	30mm
TI-102.232	32mm
TI-102.234	34mm
TI-102.236	36mm
TI-102.238	38mm
TI-102.240	40mm



2.0mm Wise-Lock Screws, Self Tapping, (Star Head)

Code	Length
TI-1131.06	6mm
TI-1131.07	7mm
TI-1131.08	8mm
TI-1131.09	9mm
TI-1131.10	10mm
TI-1131.11	11mm
TI-1131.12	12mm
TI-1131.13	13mm
TI-1131.14	14mm
TI-1131.15	15mm
TI-1131.16	16mm
TI-1131.18	18mm
TI-1131.20	20mm
TI-1131.22	22mm
TI-1131.24	24mm
TI-1131.26	26mm
TI-1131.28	28mm
TI-1131.30	30mm
TI-1131.32	32mm
TI-1131.34	34mm
TI-1131.36	36mm
TI-1131.38	38mm
TI-1131.40	40mm



- 1136-01** Self centering Double Drill Guide 2.0/1.5mm



- 1136-02** Wise-Lock Drill Sleeve 2.0mm, For Drill Bit, $\phi 1.5$ mm



- 1136-03** Bending Pin with thread For 2.0mm Wise-Lock Plates



- 1136-04** Drill Bit Quick Coupling $\phi 1.5$ mm x Length 96/82mm, 2-Flute



- 1136-05** Drill Bit Quick Coupling $\phi 2.0$ mm x Length 67/55mm, 2-Flute



- 1136-06** Countersink For 2.0/2.4/2.7mm Screws



1136-07 Screwdriver Shaft, T8 For 2.0mm Screws



1136-08 Screw Holding Sleeve For 2.0mm Cortical Screws



1136-09 Screw Holding Sleeve For 2.0mm Wise-Lock Screws



1138-09 Container with Tray and Lid For 2.0mm Implant & Instrument Set



2.4mm Wise-Lock Plate, Straight

Code	Holes	Length
TI-505.304	4	36mm
TI-505.305	5	44mm
TI-505.306	6	52mm
TI-505.307	7	60mm
TI-505.308	8	68mm
TI-505.310	10	84mm



2.4mm Adaption Plate, Straight

Code	Holes	Length
TI-821.312	12	72mm



2.4mm Wise-Lock Adaption Plate, Straight

Code	Holes	Length
TI-555.312	12	87.9mm



2.4mm Wise-Lock Condylar Plate

Code	Holes	Length
TI-559.304	Shaft 4, Head 2	37.0mm
TI-559.305	Shaft 5, Head 2	44.2mm
TI-559.306	Shaft 6, Head 2	51.4mm
TI-559.307	Shaft 7, Head 2	58.6mm



2.4mm Wise-Lock T-Adaption Plate

Code	Holes
TI-555.307T	Shaft 7, Head 2



2.4mm Wise-Lock T-Plate

Code	Holes	Length
TI-556.304	Shaft 4, Head 3	35.8mm
TI-556.305	Shaft 5, Head 3	43.0mm
TI-556.306	Shaft 6, Head 3	50.2mm
TI-556.307	Shaft 7, Head 3	57.4mm



2.4mm Wise-Lock Y-Adaption Plate

Code	Holes	Length
TI-555.304Y	Shaft 4, Head 3	38.4mm
TI-555.305Y	Shaft 5, Head 3	45.6mm
TI-555.306Y	Shaft 6, Head 3	52.8mm
TI-555.307Y	Shaft 7, Head 3	60.0mm



2.4mm Cortical Screws, Self Tapping, (Star Head)


Code	Length
TI-1144.06	6mm
TI-1144.07	7mm
TI-1144.08	8mm
TI-1144.09	9mm
TI-1144.10	10mm
TI-1144.11	11mm
TI-1144.12	12mm
TI-1144.13	13mm
TI-1144.14	14mm
TI-1144.16	16mm
TI-1144.18	18mm
TI-1144.20	20mm
TI-1144.22	22mm
TI-1144.24	24mm
TI-1144.26	26mm
TI-1144.28	28mm
TI-1144.30	30mm
TI-1144.32	32mm
TI-1144.34	34mm
TI-1144.36	36mm
TI-1144.38	38mm
TI-1144.40	40mm

2.4mm Wise-Lock Screws, Self Tapping, (Star Head)


Code	Length
TI-1145.06	6mm
TI-1145.07	7mm
TI-1145.08	8mm
TI-1145.09	9mm
TI-1145.10	10mm
TI-1145.11	11mm
TI-1145.12	12mm
TI-1145.13	13mm
TI-1145.14	14mm
TI-1145.16	16mm
TI-1145.18	18mm
TI-1145.20	20mm
TI-1145.22	22mm
TI-1145.24	24mm
TI-1145.26	26mm
TI-1145.28	28mm
TI-1145.30	30mm
TI-1145.32	32mm
TI-1145.34	34mm
TI-1145.36	36mm
TI-1145.38	38mm
TI-1145.40	40mm
TI-1145.42	42mm
TI-1145.44	44mm
TI-1145.46	46mm
TI-1145.48	48mm
TI-1145.50	50mm

1.8mm Buttress Screws, (Star Head)


Code	Length
TI-1103.112	12mm
TI-1103.114	14mm
TI-1103.116	16mm
TI-1103.118	18mm

1137-01 Self Centering Double Drill Guide 2.4/1.8mm



1137-02 Wise-Lock Drill Sleeve 2.4mm For Drill Bit Ø1.8mm



1136-06 Countersink For 2.0/2.4/2.7mm Screws



1137-03 Bending Pin With Thread For 2.4mm/2.7mm Wise-Lock Plates



1137-04 Drill Bit Quick Coupling Ø1.8mm x Length 110/82mm, 2- Flute



1137-05 Drill Bit Quick Coupling Ø2.4mm x Length 100/73mm, 2-Flute



1137-07 Screwdriver Shaft, T8 For 2.4mm/2.7mm Screws




1137-08 Holding Sleeve For 2.4mm Cortical Screws



1137-09 Holding Sleeve For 2.4/2.7mm Wise-Lock Screws



1138-10 Container with Tray and Lid For 2.4mm Implant & Instrument Set



2.7mm Wise-Lock Plate, Straight

Code	Holes	Length
TI-505.404	4	40mm
TI-505.405	5	49mm
TI-505.406	6	58mm
TI-505.407	7	67mm
TI-505.408	8	76mm
TI-505.410	10	94mm




2.7mm Wise-Lock Condylar Plate

Code	Holes	Length
TI-559.404	Shaft 4, Head 2	41.5mm
TI-559.405	Shaft 5, Head 2	49.5mm
TI-559.406	Shaft 6, Head 2	57.5mm
TI-559.407	Shaft 7, Head 2	65.5mm




2.7mm Wise-Lock T-Plate

Code	Holes	Length
TI-556.403	Shaft 3, Head 2	32.1mm
TI-556.404	Shaft 4, Head 2	39.6mm



2.7mm Wise-Lock L-Plate

Code	Holes	Length	Type
TI-813.403L	Shaft 3, Head 2	32.1mm	Left
TI-813.403R	Shaft 3, Head 2	32.1mm	Right
TI-813.404L	Shaft 4, Head 2	39.6mm	Left
TI-813.404R	Shaft 4, Head 2	39.6mm	Right



2.7mm Wise-Lock L-Plate, Oblique

Code	Holes	Length	Type
TI-1133.03L	Shaft 3, Head 2	34.6mm	Left
TI-1133.03R	Shaft 3, Head 2	34.6mm	Right
TI-1133.04L	Shaft 4, Head 2	42.1mm	Left
TI-1133.04R	Shaft 4, Head 2	42.1mm	Right



2.7mm Mini Condylar Plate

Code	Holes	Length	Type
TI-516.307L	7	57mm	Left
TI-516.307R	7	57mm	Right



2.7mm Wise-Lock Straight Plate

Code	Holes	Length
TI-1134.02	2	17mm



2.4/2.7mm Wise-Lock X-Plate

Code	Holes	Length
TI-812.301	4	22x14mm
TI-812.302	4	24x18mm
TI-812.303	4	30x20mm
TI-812.304	4	36x20mm



2.7mm Cortical Screws, Self Tapping, (Star Head)


Code	Length
TI-1242.06	6
TI-1242.08	8
TI-1242.10	10
TI-1242.12	12
TI-1242.14	14
TI-1242.16	16
TI-1242.18	18
TI-1242.20	20
TI-1242.22	22
TI-1242.24	24
TI-1242.26	26
TI-1242.28	28
TI-1242.30	30
TI-1242.32	32
TI-1242.34	34
TI-1242.36	36
TI-1242.38	38
TI-1242.40	40
TI-1242.42	42
TI-1242.44	44
TI-1242.46	46
TI-1242.48	48
TI-1242.50	50

2.7mm Wise-Lock Screws, Self Tapping, (Star Head)


Code	Length
TI-1241.06	6mm
TI-1241.07	7mm
TI-1241.08	8mm
TI-1241.09	9mm
TI-1241.10	10mm
TI-1241.12	12mm
TI-1241.14	14mm
TI-1241.16	16mm
TI-1241.18	18mm
TI-1241.20	20mm
TI-1241.22	22mm
TI-1241.24	24mm
TI-1241.26	26mm
TI-1241.28	28mm
TI-1241.30	30mm
TI-1241.32	32mm
TI-1241.34	34mm
TI-1241.36	36mm
TI-1241.38	38mm
TI-1241.40	40mm
TI-1241.42	42mm
TI-1241.44	44mm
TI-1241.46	46mm
TI-1241.48	48mm
TI-1241.50	50mm
TI-1241.52	52mm
TI-1241.54	54mm
TI-1241.56	56mm
TI-1241.58	58mm
TI-1241.60	60mm

1137-07 Screwdriver Shaft, T8 For 2.4mm/2.7mm Screws



1137-03 Bending Pin With Thread For 2.4mm/2.7mm Wise-Lock Plates



1137-09 Holding Sleeve For 2.4/2.7mm Wise-Lock Screws



1136-06 Countersink For 2.0/2.4/2.7mm Screws



1138-01 Self Centering Double Drill Guide 2.7/2.0mm



1138-02 Wise-Lock Drill Sleeve 2.7mm For Drill Bits Ø2.0mm



1138-04

Drill Bit Quick Coupling Ø2.0mm x Length 110mm, 2-Flute


1138-05

Drill Bit Quick Coupling Ø2.7mm x Length 102mm, 2-Flute


1138-41

Holding Sleeve For 2.7mm Cortical Screws


1138-11


Container with Tray and Lid For 2.7mm Implant & Instrument Set




1138.08 1.5mm Wise-Lock Hand System




1138.09 2.0mm Wise-Lock Hand System




1138.10 2.4mm Wise-Lock Hand System



1138.11 2.7mm Wise-Lock Foot System




1138.08 1.5mm Wise-Lock Hand System

Instruments For 1.5mm Wise-lock Modular Hand System

Code	Items	Qty
1135-01	Double Drill Guide 1.5/1.1mm	1
1135-02	Wise-Lock Drill Sleeve 1.5mm For Drill Bit Ø1.1mm	2
1135-03	Bending Pin With Thread For 1.5mm Wise-Lock Plates	2
1135-04	Drill Bit Quick Coupling Ø1.1mm x Length 75/61mm, 2-Flute	2
1135-05	Drill Bit Quick Coupling Ø1.5mm x length 65/50mm, 2-Flute	2
1135-06	Countersink For 1.5mm Screws	1
1135-07	Screwdriver Shaft, T4 for 1.5mm Screws	1
1135-08	Screw Holding Sleeve for 1.5mm Cortical Screws	1
1135-09	Screw Holding Sleeve for 1.5mm Wise-Lock Screws	1
1138-08	Container with Tray and Lid For 1.5mm Implant & Instrument Set	1

1138.08 1.5mm Wise-Lock Hand System
Implants For 1.5mm Wise-lock Modular Hand System


Code	Items	Qty
TI-1130.06	1.5mm Wise-Lock Screws, Self Tapping, (Star Head), Length 6mm	6
TI-1130.07	1.5mm Wise-Lock Screws, Self Tapping, (Star Head), Length 7mm	6
TI-1130.08	1.5mm Wise-Lock Screws, Self Tapping, (Star Head), Length 8mm	6
TI-1130.09	1.5mm Wise-Lock Screws, Self Tapping, (Star Head), Length 9mm	6
TI-1130.10	1.5mm Wise-Lock Screws, Self Tapping, (Star Head), Length 10mm	6
TI-1130.11	1.5mm Wise-Lock Screws, Self Tapping, (Star Head), Length 11mm	6
TI-1130.12	1.5mm Wise-Lock Screws, Self Tapping, (Star Head), Length 12mm	6
TI-1130.13	1.5mm Wise-Lock Screws, Self Tapping, (Star Head), Length 13mm	6
TI-1130.14	1.5mm Wise-Lock Screws, Self Tapping, (Star Head), Length 14mm	6
TI-1130.15	1.5mm Wise-Lock Screws, Self Tapping, (Star Head), Length 15mm	6
TI-1130.16	1.5mm Wise-Lock Screws, Self Tapping, (Star Head), Length 16mm	6
TI-1130.18	1.5mm Wise-Lock Screws, Self Tapping, (Star Head), Length 18mm	4
TI-1130.20	1.5mm Wise-Lock Screws, Self Tapping, (Star Head), Length 20mm	4
TI-1130.22	1.5mm Wise-Lock Screws, Self Tapping, (Star Head), Length 22mm	4
TI-1130.24	1.5mm Wise-Lock Screws, Self Tapping, (Star Head), Length 24mm	4
TI-1130.26	1.5mm Wise-Lock Screws, Self Tapping, (Star Head), Length 26mm	4
TI-1130.28	1.5mm Wise-Lock Screws, Self Tapping, (Star Head), Length 28mm	4
TI-1281.06	1.5mm Cortical Screws, Self Tapping, (Star Head), Length 6mm	6
TI-1281.07	1.5mm Cortical Screws, Self Tapping, (Star Head), Length 7mm	6
TI-1281.08	1.5mm Cortical Screws, Self Tapping, (Star Head), Length 8mm	6
TI-1281.09	1.5mm Cortical Screws, Self Tapping, (Star Head), Length 9mm	6
TI-1281.10	1.5mm Cortical Screws, Self Tapping, (Star Head), Length 10mm	6
TI-1281.11	1.5mm Cortical Screws, Self Tapping, (Star Head), Length 11mm	6
TI-1281.12	1.5mm Cortical Screws, Self Tapping, (Star Head), Length 12mm	6
TI-1281.13	1.5mm Cortical Screws, Self Tapping, (Star Head), Length 13mm	6
TI-1281.14	1.5mm Cortical Screws, Self Tapping, (Star Head), Length 14mm	6
TI-1281.15	1.5mm Cortical Screws, Self Tapping, (Star Head), Length 15mm	6
TI-1281.16	1.5mm Cortical Screws, Self Tapping, (Star Head), Length 16mm	6
TI-1281.18	1.5mm Cortical Screws, Self Tapping, (Star Head), Length 18mm	4
TI-1281.20	1.5mm Cortical Screws, Self Tapping, (Star Head), Length 20mm	4
TI-1281.22	1.5mm Cortical Screws, Self Tapping, (Star Head), Length 22mm	4
TI-1281.24	1.5mm Cortical Screws, Self Tapping, (Star Head), Length 24mm	4
TI-1281.26	1.5mm Cortical Screws, Self Tapping, (Star Head), Length 26mm	4
TI-1281.28	1.5mm Cortical Screws, Self Tapping, (Star Head), Length 28mm	4
TI-554.04	1.5mm Wise-Lock Plate, Straight, 4 Holes	1
TI-554.06	1.5mm Wise-Lock Plate, Straight, 6 Holes	1
TI-555.12	1.5mm Wise-Lock Adaption Plate, Straight, 12 Holes	1
TI-556.08.03	1.5mm Wise-Lock Adaption T-Plate, (3 Head Holes), 8 Holes	1
TI-556.08.04	1.5mm Wise-Lock Adaption T-Plate, (4 Head Holes), 8 Holes	1
TI-823.00R	1.5mm Wise-Lock Strut Plate, 8 Holes - Right	1
TI-823.00L	1.5mm Wise-Lock Strut Plate, 8 Holes - Left	1
TI-559.06	1.5mm Wise-Lock Condylar Plate, (2 Head Holes), 6 Holes	1
TI-559.08	1.5mm Wise-Lock Condylar Plate, (2 Head Holes), 8 Holes	1

1138.09 2.0mm Wise-Lock Hand System

Instruments For 2.0mm Wise-lock Modular Hand System

Code	Items	Qty
1136-01	Self centering Double Drill Guide 2.0/1.5mm	1
1136-02	Wise-Lock Drill Sleeve 2.0mm, For Drill Bit, Ø1.5mm	2
1136-03	Bending Pin with thread For 2.0mm Wise-Lock Plates	2
1136-04	Drill Bit Quick Coupling Ø1.5mm x Length 96/82mm, 2-Flute	2
1136-05	Drill Bit Quick Coupling Ø2.0mm x Length 67/55mm, 2-Flute	2
1136-06	Countersink For 2.0/2.4/2.7mm Screws	1
1136-07	Screwdriver Shaft, T8 For 2.0mm Screws	1
1136-08	Screw Holding Sleeve For 2.0mm Cortical Screws	1
1136-09	Screw Holding Sleeve For 2.0mm Wise-Lock Screws	1
1138-09	Container with Tray and Lid For 2.0mm Implant & Instrument Set	1

1138.09 2.0mm Wise-Lock Hand System
Implants For 2.0mm Wise-lock Modular Hand System

Code	Items	Qty
TI-102.206	2.0mm Cortical Screws, Self Tapping, (Star Head), Length 6mm	4
TI-102.207	2.0mm Cortical Screws, Self Tapping, (Star Head), Length 7mm	4
TI-102.208	2.0mm Cortical Screws, Self Tapping, (Star Head), Length 8mm	4
TI-102.209	2.0mm Cortical Screws, Self Tapping, (Star Head), Length 9mm	4
TI-102.210	2.0mm Cortical Screws, Self Tapping, (Star Head), Length 10mm	4
TI-102.211	2.0mm Cortical Screws, Self Tapping, (Star Head), Length 11mm	4
TI-102.212	2.0mm Cortical Screws, Self Tapping, (Star Head), Length 12mm	4
TI-102.213	2.0mm Cortical Screws, Self Tapping, (Star Head), Length 13mm	4
TI-102.214	2.0mm Cortical Screws, Self Tapping, (Star Head), Length 14mm	4
TI-102.216	2.0mm Cortical Screws, Self Tapping, (Star Head), Length 16mm	4
TI-102.218	2.0mm Cortical Screws, Self Tapping, (Star Head), Length 18mm	4
TI-102.220	2.0mm Cortical Screws, Self Tapping, (Star Head), Length 20mm	4
TI-102.222	2.0mm Cortical Screws, Self Tapping, (Star Head), Length 22mm	4
TI-102.224	2.0mm Cortical Screws, Self Tapping, (Star Head), Length 24mm	4
TI-102.226	2.0mm Cortical Screws, Self Tapping, (Star Head), Length 26mm	4
TI-102.228	2.0mm Cortical Screws, Self Tapping, (Star Head), Length 28mm	4
TI-102.230	2.0mm Cortical Screws, Self Tapping, (Star Head), Length 30mm	4
TI-1131.06	2.0mm Wise-Lock Screws, Self Tapping, (Star Head), Length 6mm	4
TI-1131.07	2.0mm Wise-Lock Screws, Self Tapping, (Star Head), Length 7mm	4
TI-1131.08	2.0mm Wise-Lock Screws, Self Tapping, (Star Head), Length 8mm	4
TI-1131.09	2.0mm Wise-Lock Screws, Self Tapping, (Star Head), Length 9mm	4
TI-1131.10	2.0mm Wise-Lock Screws, Self Tapping, (Star Head), Length 10mm	4
TI-1131.11	2.0mm Wise-Lock Screws, Self Tapping, (Star Head), Length 11mm	4
TI-1131.12	2.0mm Wise-Lock Screws, Self Tapping, (Star Head), Length 12mm	4
TI-1131.13	2.0mm Wise-Lock Screws, Self Tapping, (Star Head), Length 13mm	4
TI-1131.14	2.0mm Wise-Lock Screws, Self Tapping, (Star Head), Length 14mm	4
TI-1131.16	2.0mm Wise-Lock Screws, Self Tapping, (Star Head), Length 16mm	4
TI-1131.18	2.0mm Wise-Lock Screws, Self Tapping, (Star Head), Length 18mm	4
TI-1131.20	2.0mm Wise-Lock Screws, Self Tapping, (Star Head), Length 20mm	4
TI-1131.22	2.0mm Wise-Lock Screws, Self Tapping, (Star Head), Length 22mm	4
TI-1131.24	2.0mm Wise-Lock Screws, Self Tapping, (Star Head), Length 24mm	4
TI-1131.26	2.0mm Wise-Lock Screws, Self Tapping, (Star Head), Length 26mm	4
TI-1131.28	2.0mm Wise-Lock Screws, Self Tapping, (Star Head), Length 28mm	4
TI-1131.30	2.0mm Wise-Lock Screws, Self Tapping, (Star Head), Length 30mm	4
TI-505.204	2.0mm Wise-Lock Plate, Straight, 4 Holes	1
TI-505.205	2.0mm Wise-Lock Plate, Straight, 5 Holes	1
TI-505.206	2.0mm Wise-Lock Plate, Straight, 6 Holes	1
TI-505.208	2.0mm Wise-Lock Plate, Straight, 8 Holes	1
TI-555.212	2.0mm Wise-Lock Adaption Plate, Straight, 12 Holes	1
TI-555.207Y	2.0mm Wise-Lock Y-Adaption Plate, (3 Head Holes), 7 Holes	1
TI-556.207	2.0mm Wise-Lock T-Plate, (3 Head Holes), 7 Holes	1

1138.10 2.4mm Wise-Lock Hand System

Instruments For 2.4mm Wise-lock Modular Hand System

Code	Items	Qty
1137-01	Self Centering Double Drill Guide 2.4/1.8mm	1
1137-02	Wise-Lock Drill Sleeve 2.4mm For Drill Bit Ø1.8mm	2
1136-06	Countersink For 2.0/2.4/2.7mm Screws	1
1137-03	Bending Pin With Thread For 2.4mm/2.7mm Wise-Lock Plates	2
1137-04	Drill Bit Quick Coupling Ø1.8mm x Length 110/82mm, 2- Flute	2
1137-05	Drill Bit Quick Coupling Ø2.4mm x Length 100/73mm, 2-Flute	2
1137-07	Screwdriver Shaft, T8 For 2.4mm/2.7mm Screws	1
1137-08	Holding Sleeve For 2.4mm Cortical Screws	1
1137-09	Holding Sleeve For 2.4/2.7mm Wise-Lock Screws	1
1138-10	Container with Tray and Lid For 2.4mm Implant & Instrument Set	1

1138.10 2.4mm Wise-Lock Hand System
Implants For 2.4mm Wise-lock Modular Hand System

Code	Items	Qty
TI-1144.06	2.4mm Cortical Screws, Self Tapping, (Star Head), Length 6mm	4
TI-1144.07	2.4mm Cortical Screws, Self Tapping, (Star Head), Length 7mm	4
TI-1144.08	2.4mm Cortical Screws, Self Tapping, (Star Head), Length 8mm	4
TI-1144.09	2.4mm Cortical Screws, Self Tapping, (Star Head), Length 9mm	4
TI-1144.10	2.4mm Cortical Screws, Self Tapping, (Star Head), Length 10mm	4
TI-1144.11	2.4mm Cortical Screws, Self Tapping, (Star Head), Length 11mm	4
TI-1144.12	2.4mm Cortical Screws, Self Tapping, (Star Head), Length 12mm	4
TI-1144.13	2.4mm Cortical Screws, Self Tapping, (Star Head), Length 13mm	4
TI-1144.14	2.4mm Cortical Screws, Self Tapping, (Star Head), Length 14mm	4
TI-1144.16	2.4mm Cortical Screws, Self Tapping, (Star Head), Length 16mm	4
TI-1144.18	2.4mm Cortical Screws, Self Tapping, (Star Head), Length 18mm	4
TI-1144.20	2.4mm Cortical Screws, Self Tapping, (Star Head), Length 20mm	2
TI-1144.22	2.4mm Cortical Screws, Self Tapping, (Star Head), Length 22mm	2
TI-1144.24	2.4mm Cortical Screws, Self Tapping, (Star Head), Length 24mm	2
TI-1144.26	2.4mm Cortical Screws, Self Tapping, (Star Head), Length 26mm	2
TI-1144.28	2.4mm Cortical Screws, Self Tapping, (Star Head), Length 28mm	2
TI-1144.30	2.4mm Cortical Screws, Self Tapping, (Star Head), Length 30mm	2
TI-1144.32	2.4mm Cortical Screws, Self Tapping, (Star Head), Length 32mm	3
TI-1144.34	2.4mm Cortical Screws, Self Tapping, (Star Head), Length 34mm	3
TI-1144.36	2.4mm Cortical Screws, Self Tapping, (Star Head), Length 36mm	3
TI-1144.38	2.4mm Cortical Screws, Self Tapping, (Star Head), Length 38mm	3
TI-1144.40	2.4mm Cortical Screws, Self Tapping, (Star Head), Length 40mm	3
TI-1145.06	2.4mm Wise-Lock Screws, Self Tapping, (Star Head), Length 6mm	4
TI-1145.07	2.4mm Wise-Lock Screws, Self Tapping, (Star Head), Length 7mm	4
TI-1145.08	2.4mm Wise-Lock Screws, Self Tapping, (Star Head), Length 8mm	4
TI-1145.09	2.4mm Wise-Lock Screws, Self Tapping, (Star Head), Length 9mm	4
TI-1145.10	2.4mm Wise-Lock Screws, Self Tapping, (Star Head), Length 10mm	4
TI-1145.11	2.4mm Wise-Lock Screws, Self Tapping, (Star Head), Length 11mm	4
TI-1145.12	2.4mm Wise-Lock Screws, Self Tapping, (Star Head), Length 12mm	4
TI-1145.13	2.4mm Wise-Lock Screws, Self Tapping, (Star Head), Length 13mm	4
TI-1145.14	2.4mm Wise-Lock Screws, Self Tapping, (Star Head), Length 14mm	4
TI-1145.16	2.4mm Wise-Lock Screws, Self Tapping, (Star Head), Length 16mm	4
TI-1145.18	2.4mm Wise-Lock Screws, Self Tapping, (Star Head), Length 18mm	4
TI-1145.20	2.4mm Wise-Lock Screws, Self Tapping, (Star Head), Length 20mm	2
TI-1145.22	2.4mm Wise-Lock Screws, Self Tapping, (Star Head), Length 22mm	2
TI-1145.24	2.4mm Wise-Lock Screws, Self Tapping, (Star Head), Length 24mm	2
TI-1145.26	2.4mm Wise-Lock Screws, Self Tapping, (Star Head), Length 26mm	2
TI-1145.28	2.4mm Wise-Lock Screws, Self Tapping, (Star Head), Length 28mm	2
TI-1145.30	2.4mm Wise-Lock Screws, Self Tapping, (Star Head), Length 30mm	2
TI-505.304	2.4mm Wise-Lock Plate, Straight, 4 Holes	1
TI-505.306	2.4mm Wise-Lock Plate, Straight, 6 Holes	1
TI-505.308	2.4mm Wise-Lock Plate, Straight, 8 Holes	1
TI-821.312	2.4mm Adaption Plate, Straight, 12 Holes	1
TI-556.307	2.4mm Wise-Lock T-Plate, (3 Head Holes), 7 Holes	1
TI-555.307Y	2.4mm Wise-Lock Y-Adaption Plate, (3 Head Holes), 7 Holes	1

Code	Items	Qty
1137-07	Screwdriver Shaft, T8 For 2.4mm/2.7mm Screws	1
1137-03	Bending Pin With Thread For 2.4mm/2.7mm Wise-Lock Plates	2
1137-09	Holding Sleeve For 2.4/2.7mm Wise-Lock Screws	1
1136-06	Countersink For 2.0/2.4/2.7mm Screws	1
1138-01	Self Centering Double Drill Guide 2.7/2.0mm	1
1138-02	Wise-Lock Drill Sleeve 2.7mm For Drill Bits Ø2.0mm	2
1138-04	Drill Bit Quick Coupling Ø2.0mm x Length 110mm, 2-Flute	2
1138-05	Drill Bit Quick Coupling Ø2.7mm x Length 102mm, 2-Flute	2
1138-41	Holding Sleeve For 2.7mm Cortical Screws	1
1138-11	Container with Tray and Lid For 2.7mm Implant & Instrument Set	1

1138.11 2.7mm Wise-Lock Foot System
Instruments For 2.7mm Wise-lock Modular Foot System

Code	Items	Qty
TI-1242.06	2.7mm Cortical Screws, Self Tapping, (Star Head), Length 6mm	4
TI-1242.08	2.7mm Cortical Screws, Self Tapping, (Star Head), Length 8mm	4
TI-1242.10	2.7mm Cortical Screws, Self Tapping, (Star Head), Length 10mm	4
TI-1242.12	2.7mm Cortical Screws, Self Tapping, (Star Head), Length 12mm	4
TI-1242.14	2.7mm Cortical Screws, Self Tapping, (Star Head), Length 14mm	4
TI-1242.16	2.7mm Cortical Screws, Self Tapping, (Star Head), Length 16mm	4
TI-1242.18	2.7mm Cortical Screws, Self Tapping, (Star Head), Length 18mm	4
TI-1242.20	2.7mm Cortical Screws, Self Tapping, (Star Head), Length 20mm	4
TI-1242.22	2.7mm Cortical Screws, Self Tapping, (Star Head), Length 22mm	1
TI-1242.24	2.7mm Cortical Screws, Self Tapping, (Star Head), Length 24mm	1
TI-1242.26	2.7mm Cortical Screws, Self Tapping, (Star Head), Length 26mm	1
TI-1242.28	2.7mm Cortical Screws, Self Tapping, (Star Head), Length 28mm	1
TI-1242.30	2.7mm Cortical Screws, Self Tapping, (Star Head), Length 30mm	1
TI-1242.32	2.7mm Cortical Screws, Self Tapping, (Star Head), Length 32mm	3
TI-1242.34	2.7mm Cortical Screws, Self Tapping, (Star Head), Length 34mm	3
TI-1242.36	2.7mm Cortical Screws, Self Tapping, (Star Head), Length 36mm	3
TI-1242.38	2.7mm Cortical Screws, Self Tapping, (Star Head), Length 38mm	3
TI-1242.40	2.7mm Cortical Screws, Self Tapping, (Star Head), Length 40mm	3
TI-1241.06	2.7mm Wise-Lock Screws, Self Tapping, (Star Head), Length 6mm	4
TI-1241.08	2.7mm Wise-Lock Screws, Self Tapping, (Star Head), Length 8mm	4
TI-1241.10	2.7mm Wise-Lock Screws, Self Tapping, (Star Head), Length 10mm	4
TI-1241.12	2.7mm Wise-Lock Screws, Self Tapping, (Star Head), Length 12mm	4
TI-1241.14	2.7mm Wise-Lock Screws, Self Tapping, (Star Head), Length 14mm	4
TI-1241.16	2.7mm Wise-Lock Screws, Self Tapping, (Star Head), Length 16mm	4
TI-1241.18	2.7mm Wise-Lock Screws, Self Tapping, (Star Head), Length 18mm	4
TI-1241.20	2.7mm Wise-Lock Screws, Self Tapping, (Star Head), Length 20mm	4
TI-1241.22	2.7mm Wise-Lock Screws, Self Tapping, (Star Head), Length 22mm	2
TI-1241.24	2.7mm Wise-Lock Screws, Self Tapping, (Star Head), Length 24mm	2
TI-1241.26	2.7mm Wise-Lock Screws, Self Tapping, (Star Head), Length 26mm	2
TI-1241.28	2.7mm Wise-Lock Screws, Self Tapping, (Star Head), Length 28mm	2
TI-1241.30	2.7mm Wise-Lock Screws, Self Tapping, (Star Head), Length 30mm	2
TI-505.404	2.7mm Wise-Lock Plate, Straight, 4 Holes	1
TI-505.406	2.7mm Wise-Lock Plate, Straight, 6 Holes	1
TI-556.404	2.7mm Wise-Lock T-Plate, (2 Head Holes), 4 Holes	1
TI-813.403L	2.7mm Wise-Lock L-Plate, (2 Head Holes), 3 Holes - Left	1
TI-813.403R	2.7mm Wise-Lock L-Plate, (2 Head Holes), 3 Holes - Right	1
TI-813.404L	2.7mm Wise-Lock L-Plate, (2 Head Holes), 4 Holes - Left	1
TI-813.404R	2.7mm Wise-Lock L-Plate, (2 Head Holes), 4 Holes - Right	1
TI-559.407	2.7mm Wise-Lock Condylar Plate, (2 Head Holes), 7 Holes	1

1138-12

Handle with Mini Quick Coupling


1138-13

Handle With Quick Coupling (Straight)


1138-14

Depth Gauge measuring upto 30mm


1138-15

Depth Gauge measuring upto 50mm


1138-17

Screw Holding Forceps


1138-18

Cutting Forceps



1138-23

Torque Limiting attachment, 0.8Nm, with AO quick coupling (Optional)


1138-24

Torque Limiting attachment, 0.4Nm, with AO quick coupling (Optional)


1138-26

Hohmann Retractor, Width 15mm


1138-27

Hohmann Retractor, Width 6.5mm


1138-28

Periosteal Elevator, Round


1138-29

Periosteal Elevator, Flat



1138-30

Cutting Plier (For pins and plates)


1138-31

Bending Pliers


1138-33

Reduction Forceps, Pointed, Small


1138-36

Plate File


1138-37

Reduction Forceps, Serrated Jaws


1138-38

Star Screw Driver for 1.5mm Screws



1138-39


Star Screw Driver for 2.0mm Screws



1138-40

Star Screw Driver For 2.4/2.7mm Screws


1138-35

Container for Hand & Foot Instrument & Implant set



5001.000 Instrument System for Modular Hand & Foot System


Code	Items	Qty
1138-08	Container with Tray and Lid For 1.5mm Implant & Instrument Set	1
1138-09	Container with Tray and Lid For 2.0mm Implant & Instrument Set	1
1138-10	Container with Tray and Lid For 2.0mm Implant & Instrument Set	1
1138-11	Container with Tray and Lid For 2.7mm Implant & Instrument Set	1
1138-12	Handle with Mini Quick Coupling	1
1138-13	Handle With Quick Coupling (Straight)	1
1138-14	Depth Gauge measuring upto 30mm	1
1138-15	Depth Gauge measuring upto 50mm	1
1138-17	Screw Holding Forceps	1
1138-18	Cutting Forceps	1
1138-23	Torque Limiting attachment, 0.8Nm, with AO quick coupling (Optional)	1
1138-24	Torque Limiting attachment, 0.4Nm, with AO quick coupling (Optional)	1
1138-26	Hohmann Retractor, Width 15mm	1
1138-27	Hohmann Retractor, Width 6.5mm	1
1138-28	Periosteal Elevator, Round	1
1138-29	Periosteal Elevator, Flat	1
1138-30	Cutting Plier (For pins and plates)	1
1138-31	Bending Pliers	2
1138-33	Reduction Forceps, Pointed, Small	1
1138-36	Plate File	1
1138-37	Reduction Forceps, Serrated Jaws	1
1138-38	Star Screw Driver for 1.5mm Screws	1
1138-39	Star Screw Driver for 2.0mm Screws	1
1138-40	Star Screw Driver For 2.4/2.7mm Screws	1
1138-35	Container for Hand & Foot Instrument & Implant set	1



USA

Auxein Inc.

1500 Nw 89th Court, Suite 205, Doral, Florida 33172

Tel: +1 305 395 6062

E Fax: +1 305 395 6262

Email: USoffice@auxein.com

INDIA

Auxein Medical Pvt. Ltd.

Plot No. 168-169-170, Phase-4, Kundli Industrial Area,
HSIIDC, Sector-57, Sonapat - 131028, Haryana

Tel: +91 99106 43638

Fax: +91 86077 70197

Email: info@auxein.com

www.auxein.com