



ic-cerclage



implantcast


ic-cerclage

**material : pure titanium acc. to
ISO 5832-2**

REF 0060-1008 Ø 8mm
REF 0060-1018 Ø 8mm



**ic-cerclage container
REF 7999 0060**



**ic-cerclage tensioner
REF 0060 1010**



**ic-strap guide 10°
REF 0060 1020**



**ic-strap guide 80°
REF 0060 1021**



**ic-punch
REF 0060 1030**



ic-titanium cerclage wire



**material : pure titanium acc. to
ISO 5832-2**

REF 0060-2014 Ø 1,4 mm, length 5m
REF 0060-2013 Ø 1,4 mm, length 50cm



implantcast GmbH
Lüneburger Schanze 26
D-21614 Buxtehude
Germany
phone: +49 4161 744-0
fax: +49 4161 744-200
e-mail: info@implantcast.de
internet: www.implantcast.de



Your local distributor:

